




ABOUT BREATHE-ON DEVICE



-  (KC certification) - Korea national integrated certification for safety, health, environment and quality
-  (FCC certification) - US Federal Communications Commission
-  (CE certification) - European Union (EU) integrated standards for safety, health, environment and consumer protection
-  (CB certification) - International certification for safety of electric and electronic products and electromagnetic waves(required for international trades)
-  (Selected by Good Design 2017) - Korea Design Promotion Agency
- Aptness assessment - Korea University Hospital
- Usability assessment - Senior Experience Center at Kyemyung University, Daegu
- Functionality Test Result - Daegu Mechatronics & Materials Institute
- Ergonomic Design Awards - Korean Agency for Technology and Standards

ABOUT BREATHE-ON COMPONENTS



01. Measurement / evaluation program

» Measuring and evaluating the user's breathing condition and breathing ability help proper breathing exercise.

02. Customized training program

» Different difficulty levels for the user's breathing condition and breathing ability

- » Data record of inspiration, expiration, time, and strength
- » Various breathing patterns by 12 different contents
- » Connecting with a microphone instead of the breathing device, you can play vocalizations.



JOIN US WITH THE BREATHE-ON

Breathing / Vocalization / Cognition development / Learning / Game




Excellent product
registered with PPS

Eco-friendly
product package
Safe product design
KC Certified

Oleophobic
(anti smudge) coating
on the glass a low
coefficient of friction
4mm-thick
tempered glass

32-inch touch



Safe product design
to prevent scrapes on
the foot- the legs of
the product located
underneath the screen
table

Rounded
corners



Various contents

Various multi-touch contents
that children can enjoy and study
altogether




10 points

Simultaneously available for
6 users



Android

Simultaneous multi-touch
function in Android environment



- » 12 contents for breathing, vocal training and exercise
- » Special contents to enhance breathing ability and lung capacity.



- » Very special contents of more than 50 varieties that people of all ages can enjoy learning, games, drawing, playing musical instrument, memorizing, maze and puzzle, using visual, auditory and tactile sense

ABOUT THE DISPLAY AND CONVENIENT DEVICE



Smartphone
(Android)



Tablet
(Android)



Touch monitor
(Android)



Mobile stand



UV sterilizer

* Currently only available for Android (iOS is under development)

* The above images for reference only-the actual product image can be different.




APPLICATIONS



Constant increase in users


Big-Data Generation

Those consent with information use
|
Providing health information



EXPECTED EFFECTS AND BUSINESS DIRECTION

- Sustainable Product Development
- New products compatible with the related products
- Various content development




JOIN US WITH THE BREATHE-ON

Breathing / Vocalization / Cognition development / Art / Game

Special features of the Feeling Pen

- » Use anytime and anywhere with your mobile phone or tablet PC
- » Easy to catch attention of children more familiar with the smartphone than paper.
- » Vibration feedback system for children to correct their drawing while enjoying tactile sensation.
- » Interactive learning with various educational contents.
- » Free drawing
- » Brain stimulation learning method with vibration!



4D Haptic Brain



Designing



Drawing



Delighting



Dreaming



A variety of designs!

You can choose any petite character you like.
Show off and change at any time!



JOIN US WITH THE BREATHE-ON

Breathing / Vocalization / Cognition development / Art / Game

Domestic technology

Vibration feedback

Vibration feedback technology dramatically changes the paradigm of a digital pen.

A new sensation for children learning, seniors worrying about dementia and those with hearing and visual impairments



Vibration Feedback Point

Drawing line

Guideline

Vibration Feedback Point

Tactile language using vibration



JOIN US WITH THE BREATHE-ON

Breathing / Vocalization / Cognition development / Art / Game

01 Seven-phase art

<Art with Oringa> Art education contents lead students to learn naturally a sense of color using 24 colors and systematic drawing techniques.



Familiarize yourself with the names of 24 colors via color matching.



Learn the color combination while coloring.



Learn how to draw from drawing along the line



Learn how to express feelings of an object.



Learn how to draw a movement of an object.



Learn how to draw using numbers, alphabets, and consonants.



Express your feelings freely while sketching.



Enjoy Oringa animation.

JOIN US WITH THE BREATHE-ON

Breathing / Vocalization / Cognition development / Art / Game


02 Measurement function for learning activity

Identify favorite colors of the children and use frequency of them in drawing.



04 Creative education

Increasing creativity and self-esteem via the Dream Sketch!



03 Multilingual support

Supporting 6 languages(Korean, English, Chinese, Japanese, Spanish and Russian)



05 Printing system

You can print out the contents and draw on the paper.



JOIN US WITH THE BREATHE-ON

Breathing / Vocalization / Recovery from fatigue and pain on muscle and joint

Tactile sensor

When the main unit is taken off, the stimulation is immediately stopped and the alarm system is activated for safety.

*Any time anywhere!
Easy to carry!*



At home



For the break time



Watching TV



Reading



Before and after exercise



In the office

Electrical discharge to remove residual current

Residual current is discharged at power off to prevent the electric shock.

- Easy to use
- Four stimulation modes
Toning / Firming / Bulking / Exercise
- Easily Attachable and Detachable
- Built-in lithium polymer rechargeable battery (Needs to be charged for use)
- Excellent portability (special case)

Soft start

Even if the stimulation is interrupted, the main unit memorizes the stimulation mode and setting level. When reconnecting, it goes up from the level 0 to the existing setting level (1 level per second).



* This product is a medical device.

JOIN US WITH THE BREATHE-ON

Breathing / Vocalization / Recovery from fatigue and pain on muscle and joint



* This product is a medical device.

Pad and gel (hot and cool) for waist, knee, and elbow- separately purchasable